

**Figure 1.** Thermalytix system generated images for case ID B146 showing malignant breast lesions. Images were captured from five different angles: frontal (A), lateral (Right-B, Left-D), and oblique (Right-C, Left- E) views. The regions with pink borders denote warm spots whereas blue borders denote suspicious malignant lesions.

Using the hotspot, areolar, vascular, and overall thermal and/or demographic properties, Thermalytix generated thermobiological, areolar, vascular, and ensemble scores of 0.32, 0.02, 0.69, and 0.66, respectively. A score >0.5 indicated a high likelihood of abnormal lesions. These scores are further combined to obtain a B-Score of 4 indicating a high likelihood of a cancerous lesion in the right breast and a breast ultrasound was recommended for further diagnosis. Table 1 indicates the significance of the B Score and a comparison with

**Table 1.** B-Score description.

B-SCORE	THERMALYTIX INTERPRETATION	ACR BI-RADS ASSESSMENT CLASSIFICATION
B-Score 0	Inconclusive. Repeat Thermalytix test, probably an error in thermal image capture.	Category 0: Incomplete – need additional imaging evaluation and/or prior mammograms for comparison
B-Score 1	Normal. No evident signs of breast malignancy in the thermal image.	Category 1: Negative
B-Score 2	Normal. No evident signs of breast malignancy in the thermal image. However, some regions of thermal activity are observed.	Category 2: Benign
B-Score 3	Borderline chances of malignancy. Suggest repeating the Thermalytix test in 3-6 months. For patients with a family cancer history and other risk factors, suggest a diagnostic workup.	Category 3: Probably benign
B-Score 4	Suspicious for malignancy, suggest diagnostic workup	Category 4: Suspicious • Category 4A: Low suspicion for malignancy • Category 4B: Moderate suspicion for malignancy • Category 4C: High suspicion for malignancy
B-Score 5	Highly suspicious for malignancy, suggest diagnostic workup	Category 5: Highly suggestive of malignancy


the American College of Radiology-Breast Imaging Reporting and Data Systems (ACR BI-RADS) assessment category.

The patient underwent a bilateral mammogram and USG of the right breast lesion. Bilateral mammography evaluation categorized breast density as BI-RADS A parenchymal pattern and detected a well-defined lesion measuring “20 × 20 mm” in the upper outer quadrant of the right breast at 10 O’clock position with mildly spiculated margins and with few microcalcifications. The right axillary nodes were enlarged, with the largest node measuring “18 × 12 mm”. The final impression of the lesion was categorized as ACR BI-RADS 4. USG of the right breast revealed spiculated hypoechoic mass lesion measuring “16 × 12 mm” in the 12 O’clock position Circle II Zone B of malignant etiology and was categorized as BI-RADS 4C lesion - moderate to a substantial likelihood of malignancy 51%-95%.

The histopathology of a wide local excision specimen of the right breast mass showed invasive carcinoma - no special type with lymphovascular invasion. The tumor was grade 2, measuring “20 × 18 × 18 mm” and the resection margins were clear. Dissection of the axillary lymph nodes showed micrometastasis in one of five lymph nodes. The pathological staging for this patient was pT1cN1mi. Immunohistochemical evaluation revealed estrogen receptor positivity, progesterone receptor positivity, HER 2 negativity, and Ki-67 proliferative index up to 25%.

**Case 2:** A 46-year-old premenopausal woman presented to a tertiary care hospital with a lump in the upper inner quadrant (UIQ) of her right breast for 2 months. She was the mother of four breastfed children and had no family history of cancer.

The patient underwent the Thermalytix test. Three hotspots were identified that extended to 45.9% of the UIQ of the right breast. The hotspots were irregular and distorted with a 1.14°C increase with respect to the surrounding



**Figure 2.** Thermalytix system generated images for case ID: M33 showing malignant breast lesions. Images were captured from five different angles: frontal (A), lateral (Right-B, Left-D), and oblique (Right-C, Left-E) views. The regions outlined with pink borders denote warm spots whereas regions outlined with blue borders denote suspicious malignant lesions.

region. There was neither thermal nor areolar symmetry (0%) and six major vessels with five vessel branches are seen on the right breast.

The high Thermalytix risk scores (TRSs; thermobiological: 0.85, areolar: 0.95, vascular: 0.63, and ensemble: 0.91) and a high B-score of five indicated a high suspicion of a cancerous lesion in the right breast. The AI tool also generated annotated thermal images that showed multiple focal thermal patterns with a thermal increase near the areolar region in the right breast (Figure 2).

Bilateral mammographic evaluation revealed that both breasts were composed of mixed fat and moderate fibro glandular parenchyma. The breast density was categorized as BI-RADS type C. There was no evidence of focal lesion, micro- or macro-calcification, skin thickening, or nipple retraction in both right and left breasts. Mammographic impression was provided as ACR BI-RADS 0, recommending further investigation. An ultrasound guided

fine needle aspiration cytology of the right breast lump revealed paucicellular smears. A core biopsy of the right breast was performed to further investigate the lump. The pathology examination results reported the lesion as grade 2 IDC without angioinvasion.

Thus, this was a case of malignancy that was accurately identified by Thermalytix although it was BI-RADS 0 on mammography.

These are unique cases where a malignant breast lesion was detected by Thermalytix and was confirmed on histopathological examination. Table 2 summarizes the characteristics of two cases and Table 3 summarizes the thermal characteristics of the hot lump presented in this report. Thus, a comprehensive Thermalytix-based approach has the potential to identify malignant breast lesions such as IDC.

**Discussion**

Mammography is the WHO-recommended modality for screening in developed countries [5]. However, mammography may not be the most appropriate screening modality in low- and middle-income countries (LMIC). Mammography also has low sensitivity in women with dense breasts and is usually complemented by USG or magnetic resonance imaging to increase the early identification of malignancies in dense breasts [6]. As these tests are resource-intensive and expensive, the benefits and risks of screening mammography have been debated in LMICs.

Thermography is a medical imaging technique that records changes in the surface temperature of a human body using infrared radiation generated by the body’s surface [7]. Currently available diagnostic modalities (mammography and USG) only detect morphologic changes in the breast once it has grown to a particular size, whereas thermography identifies physiological changes in the breast tissue, such as temperature and blood flow, that are indicative of tumor growth [7]. Thermography, a non-invasive method of breast imaging, was approved as an adjunct imaging modality to mammography in 1982 by The US Food and Drug Administration. However, thermography was not widely accepted due to the complexity involved in interpreting thermograms with the naked eye and resulting in a high false-positivity rate [8]. In recent years, infrared breast thermography is re-emerging as a modality with promising results for detecting breast

**Table 2.** Characteristics of clinical cases.

CASE ID	AGE	SYMPTOMS	MENO-PAUSE	BREAST DENSITY	THERMALYTX	MAMMOGRAPHY	HISTOPATHOLOGY
B146	47	Lump	Yes	A	B-Score 4positive	BI-RADS 4 (Suspicious of malignancy)	IDC Grade 2
M33	46	Lump	No	C	B-Score 5 positive	BI-RADS 0 (Inconclusive)	IDC Grade 2

**Table 3.** Thermal characteristics of the hot lump.

THERMALYTIX SCORES	B146	M33
Thermo biological score <sup>a</sup>	0.32	0.85
Areolar score <sup>a</sup>	0.02	0.95
Vascular score <sup>a</sup>	0.69	0.63
Ensemble score <sup>a</sup>	0.66	0.91
Temperature	1.32°C increase wrt surrounding region- hot lump	1.14°C increase wrt surrounding region- hot lump

<sup>a</sup>score >0.5 indicates a high likelihood of the presence of abnormal lesions.

**Table 4.** Comparison of multiple AI-based breast screening solutions.

MEDICAL SOFTWARE DEVICE	PERFORMANCE	REFERENCE
Thermalytix Niramai health analytix	Sensitivity 91.02% Specificity 82.39%	Kakileti et al. [10]
Transpara™ (ScreenPoint Medical/ Incepto™)	Sensitivity 86% Specificity 79%	Rodríguez-Ruiz et al. [11]
Mammography intelligent assessment®	Sensitivity 85%-97% Specificity 50%-94%	Sharma et al. [12].
Therapixel Nice, France	Sensitivity 69% Specificity 73%	Pacilè et al. [13]

cancer because of better thermal cameras, improved algorithms, and reduced manual error by using AI [9].

In our technique, 31 features are extracted from the hot spot region of interest. And to differentiate the malignant heat patterns, we use a random forest (RF) classifier with 200 decision trees. To assess the vascular score, we extract 21 features from vessels detected in both the breast regions. These features are then fed to an RF classifier with 200 decision trees to calculate the Ensemble score. To obtain the final TRS, we combine the proposed vascular score, and hotspot score with critical symptoms like lump and nipple discharge each with different weightages which is explained in the study of Kakileti et al. [4].

Table 4 depicts the comparison of different AI-based screening solutions in terms of their performance. In recent times, thermography has been recommended as a low-cost imaging modality and an alternative to mammography for the early detection of breast cancer lesions in economically backward rural settings [14].

Thermalytix, an AI-based thermography screening tool, was compared to standard of care diagnostics in a multi-site clinical trial of women with symptomatic breast lesions and found to be non-inferior to mammography. Multiple studies have shown the overall sensitivity and

specificity of Thermalytix are high [3,10]. The above-described cases depict the technical advances in thermography, coupled with machine learning to improve its diagnostic ability and use as a potential adjunctive tool along with standard imaging modalities.

## Conclusion

IDC is the most frequent type of breast cancer and has a high survival rate when diagnosed and treated early. Thermalytix is a portable, radiation-free breast cancer screening technology that uses AI to analyze thermal images. Clinical trials have shown that it is non-inferior to mammography and is a potential screening solution in LMICs. Thermalytix is an automated, affordable, privacy-aware solution for women of all age groups. The analysis of the above two clinical cases of IDC suggests that Thermalytix has the potential in detecting invasive breast cancer in women undergoing breast screening.

## What is new?

The presented work evaluates a new technique for the early detection of breast cancer, called Thermalytix, which is a computer-aided diagnostic solution that uses artificial intelligence algorithms over thermal images to detect possible malignant lesions. In this work, the authors describe two clinical cases of symptomatic female patients with invasive ductal carcinoma accurately detected by Thermalytix and further evaluated and confirmed by standard of care imaging and tissue diagnosis methods. The authors believe that this case report will be of interest to your readers because Thermalytix is a novel screening tool that is low-cost, accurate, automated, and portable, allowing it to be used in any clinic.

## List of Abbreviations

AI	Artificial intelligence
BI-RADS	Breast imaging reporting & data system
LMICs	Low- and middle-income countries
MRI	Magnetic resonance imaging
USG	Ultrasonography
UIQ	Upper inner quadrant

## Conflict of interest

C.G., L.K., S.T.K., H.M., S.C., and G.M. are employees of Niramai Health Analytix Private Ltd and have direct conflict of interest. A.S., V.B., R.B., and B.A. do not have any conflict of interest regarding the publication of this article.

## Funding

None.

## Consent for publication

Informed written consent was taken from the patients to participate in the study.

## Ethical approval

Not required.

### Author details

Charitha Gangadharan<sup>1</sup>, Lakshmi Krishnan<sup>1</sup>, Siva Teja Kakileti<sup>1</sup>, Sathiakar Collison<sup>2</sup>, Akshita Singh<sup>3</sup>, Venkatraman Bhat<sup>3</sup>, Richa Bansal<sup>6</sup>, Bharat Aggarwal<sup>6</sup>, Geetha Manjunath<sup>1</sup>

1. Niramai Health Analytix Pvt. Ltd., Bangalore, India
2. Sriram Cardio-Neuro Center, Kalra Hospital, New Delhi, India
3. Mazumdar Shaw Medical Centre, Narayana Hrudayalaya Health City, Bangalore, Karnataka, India
4. Narayana Hrudayalaya Limited, Bangalore, India
5. Narayana Hrudayalaya Health City, Bangalore, India
6. Max Super Speciality Hospital, New Delhi, India

### References

1. Makki J. Diversity of breast carcinoma: histological subtypes and clinical relevance. *Clin Med Insights Pathol*. 2015;8:23–31. <https://doi.org/10.4137/CPATH.S31563>
2. Parambil NA, Philip S, Tripathy JP, Philip PM, Duraisamy K, Balasubramanian S. Community engaged breast cancer screening program in Kannur District, Kerala, India: a ray of hope for early diagnosis and treatment. *Indian J Cancer*. 2019;56(3):222–7. [https://doi.org/10.4103/ijc.IJC\\_397\\_18](https://doi.org/10.4103/ijc.IJC_397_18)
3. Singh A, Bhat V, Sudhakar S, Namachivayam A, Gangadharan C, Pulchan C, et al. Multicentric study to evaluate the effectiveness of Thermalytix as compared with standard screening modalities in subjects who show possible symptoms of suspected breast cancer. *BMJ Open*. 2021;11(10):e052098. <https://doi.org/10.1136/bmjopen-2021-052098>
4. Kakileti ST, Madhu HJ, Manjunath G, Wee L, Dekker A, Sampangi S. Personalized risk prediction for breast cancer pre-screening using artificial intelligence and thermal radiomics. *Artif Intell Med*. 2020;105:101854. <https://doi.org/10.1016/j.artmed.2020.101854>
5. WHO. Position paper on mammography screening. Geneva, Switzerland: World Health Organization; 2014.
6. Mann RM, Athanasiou A, Baltzer PAT, Camps-Herrero J, Clauser P, Fallenberg EM, et al. European society of breast imaging (EUSOBI). breast cancer screening in women with extremely dense breasts recommendations of the European society of breast imaging (EUSOBI). *Eur Radiol*. 2022;32(6):4036–45. <https://doi.org/10.1007/s00330-022-08617-6>
7. Lahiri BB, Bagavathiappan S, Jayakumar T, Philip J. Medical applications of infrared thermography: a review. *Infrared Phys Technol*. 2012;55(4):221–35. <https://doi.org/10.1016/j.infrared.2012.03.007>
8. Malik AM. *New Perspectives in Breast Imaging*. London: IntechOpen. 2017:184. Available from: <https://www.intechopen.com/books/6047> doi: 10.5772/67606
9. Ghobadi H, Thainimit S, Sugino N, Gansawat D, Zadeh HG. Comparative accuracy of digital infra-red thermal imaging (DITI) in breast cancer diagnosing. *J Chem Pharm Res*. 2016;8:577–83.
10. Kakileti ST, Madhu HJ, Krishnan L, Manjunath G, Sampangi S, Ramprakash HV. Observational study to evaluate the clinical efficacy of Thermalytix for detecting breast cancer in symptomatic and asymptomatic women. *JCO Glob Oncology*. 2020;6:1472–80. <https://doi.org/10.1200/GO.20.00168>
11. Rodríguez-Ruiz A, Krupinski E, Mordang JJ, Schilling K, Heywang-Köbrunner SH, Sechopoulos I, et al. Detection of breast cancer with mammography: effect of an artificial intelligence support system. *Radiology*. 2019;290(2):305–14. <https://doi.org/10.1148/radiol.2018181371>
12. Sharma N, Ng AY, James JJ, Khara G, Ambrozay E, Austin CC, et al. Large-scale evaluation of an AI system as an independent reader for double reading in breast cancer screening. *medRxiv*. 2021. <https://doi.org/10.1101/2021.02.26.21252537>
13. Pacilè S, Lopez J, Chone P, Bertinotti T, Grouin JM, Fillard P. Improving breast cancer detection accuracy of mammography with the concurrent use of an artificial intelligence tool. *Radiol Artif Intell*. 2020;2(6):e190208. <https://doi.org/10.1148/ryai.2020190208>
14. Khan AA, Arora AS. Thermography as an economical alternative modality to mammography for early detection of breast cancer. *J Healthc Eng*. 2021;2021. <https://doi.org/10.1155/2021/5543101>

### Summary of the case

1	<b>Patient (Gender and age)</b>	Female, 47 years old	Female, 46 years old
2	<b>Final diagnosis</b>	IDC-Grade 2	IDC-Grade 2
3	<b>Symptoms</b>	Yes	Yes
4	<b>Medications</b>	None	None
5	<b>Clinical procedures</b>	Wide excision	Core biopsy
6	<b>Specialty</b>	Radiology and imaging	Radiology and imaging